

Technical Data Sheet

P-DERM®

PS-1862

Product Description:


P-DERM® PS-1862 is a trilaminate consisting of high adhesion silicone gel adhesive, polyurethane film and medical grade acrylic adhesive where consistent adhesion and repositionability is required. PS-1862 is supplied with a 100 micron polyethylene liner on the silicone gel adhesive side and 54# bleached kraft siliconized liner on the medical grade adhesive side of the construction.

Construction / Properties:

Property	Value	Test Method
Thickness (excluding liners)	0.22mm	QSP-726
Silicone Gel Coat Weight	150 gsm	QSP-724
Carrier Type	Polyurethane Film	
Carrier Thickness	25 micron	QSP-726
Adhesive Coat Weight	37 gsm	QSP-724
Silicone Adhesion	2.5 N/25mm	QSP-723
Medical Grade Adhesive Adhesion	>13 N/25mm	QSP-722
Continuous Use Conditions	10 - 70 C	QSP-754
Biocompatibility	Cytotoxicity: Pass Primary Dermal: Pass Buehler Sensitization: Pass	ISO 10993
Recommended Sterilization Method	Ethylene Oxide	

Features:

- Excellent Instantaneous Tack
- Atraumatic Removal from Skin
- Hypoallergenic
- · Ability to be Removed and Repositioned
- Excellent Converting Properties



Specific tests should be performed by the end user to determine the product stability for the particular application.

For Additional Information:

E-mail: sales@polymerscience.com

Toll Free: +I 888.533.7004

Web: www.polymerscience.com

Revision: 013124

