

# Technical Data Sheet

P-DERM®

# PS-1818

### **Product Description:**


P-DERM® PS-1818 is a silicone gel adhesive tape consisting of a high adhesion silicone gel skin contact adhesive coated onto a polyurethane carrier. PS-1818 is supplied with a 75 micron matte polycarbonate.

## **Construction / Properties:**

Property	Value	Test Method
Total Thickness	0.27mm	QSP-726
Silicone Gel Coat Weight	250 gsm	QSP-724
Polyurethane Carrier	28 micron	QSP-724
Silicone Gel Adhesion	1.8 N/25mm	QSP-723
Biocompatibility	Cytotoxicity: Pass Primary Dermal: Pass Buehler Sensitization: Pass	ISO 10993
Continuous Use Range	10 - 70 C	WI-8.2-13
Recommended Sterilization Method	Ethylene Oxide	

#### **Features:**

- Excellent Instantaneous Tack
- Atraumatic Removal from Skin
- Hypoallergenic
- Ability to be Removed and Repositioned
- Excellent Converting Properties



Specific tests should be performed by the end user to determine the product stability for the particular application.

#### For Additional Information:

E-mail: sales@polymerscience.com

Toll Free: +I 888.533.7004

Web: www.polymerscience.com

Revision: 062718

